

Line ECO

Laboratory drying ovens and incubators

NEW



MMM Group



- New LCD panel 9 inch
- resistant new foil keypad with pleasant Softtouch surface
- clear mechanical keypad response
- transfective FSTN display with COG technology (it is highlighted + it uses the external lighting reflexion – higher intensity of external lighting improves the display readability)
- brightened keys and logo integrated in the foil keypad
- adjustable display contract (e.g. in case of the device to be placed too high or too low)
- large long distance visible display
- main chamber data (T/Rh) on the screen
- intuitive control of the device
- integrated datalogger as option
- three wires PT 100 (protective + control)
- patented userfriendly door lock
- ergonomic door with reliable double closure
- power supply 115/230 V
- Fuzzy logic control system
- vertical air flow in the chamber
- vertical chamber construction saving space in the laboratory
- stainless steel chamber for long term stable chamber conditions
- double jacket chamber for exact temperature regulation
- easy cleaning thanks to simply removable panels and rounded chamber corners



Available from 2017



MMM Medcenter Einrichtungen GmbH
Sammelweisstrasse 6
D-82152 Planegg / München
Germany

tel: +49 89 899 22 620
fax: +49 89 899 22 630
medcenter@mmm-medcenter.de
www.mmm-medcenter.de



protecting human health

ECO Flyer_052016_EN/PR